



### TESTIMONIAL

"Rodents will attack the bait as efficiently with a plastic wrap and the bait LASTS A LOT LONGER in the wrapper. From a cost analysis perspective for the pest controller, replacing baits regularly due to other pests costs money, but in this regard we find other pests go away and we only get hits from the intended rodent target."

**Charles Said-Vassalo**  
Commercial Pest Management



Liphatech


Product information contact:  
Jerome Beline +61 421 946 857  
Romain Broch +61 416 191 594  
email: [info@liphatech.com.au](mailto:info@liphatech.com.au)

Disclaimer: Liphatech S.A.S. shall not be liable for any consequential or other loss or damage relating to the supply or subsequent handling or use of this product, unless such liability by law cannot be lawfully excluded or limited. All warranties, conditions or rights implied by statute or other law which may be lawfully excluded are so excluded. Where the liability of Liphatech S.A.S. for breach of any such statutory warranties and conditions cannot be lawfully excluded but may be limited to it re-supplying the product or an equivalent product or the cost of a product or an equivalent product, then the liability of Liphatech S.A.S. for any breach of such statutory warranty or condition is so limited.

® Registered trademark of Liphatech S.A.S

ALWAYS READ THE PRODUCT LABEL BEFORE USE

Available from:



# MAKI® BLOCK

Available as



**MAKI® ORIGINAL**  
IF YOU PREFER THEM **UNDRESSED**  
(NAKED)

or



**MAKI® WRAPPED**  
WE'VE GOT YOU COVERED

### PRACTICALITY FEATURES

- MAKI® Block has a central hole so it can be easily secured inside a bait station
- High number of baits per bucket (375 for MAKI® Original, 275 for MAKI® Wrapped) allows more baiting points at times of low activity
- Convenient rectangular bucket for optimum storage

### SAFETY FEATURES

- Bromadiolone is an anti-coagulant. Antidote Vitamin K1 is effective in case of accidental poisoning of children or pets
- MAKI Block contains Bitrex, a bittering agent that reduces the risk of consumption by children or pets



MAKI® Original  
7.5kg - 375 Blocks

MAKI® Wrapped  
5.5kg - 275 Blocks

**LIPHA**TECH  
[www.liphatech.com.au](http://www.liphatech.com.au)

**LIPHA**TECH  
[www.liphatech.com.au](http://www.liphatech.com.au)



# MAKI ORIGINAL

IF YOU PREFER THEM UNDRRESSED  
(NAKED)


# MAKI WRAPPED

WE'VE GOT YOU COVERED



## YOUR EVERYDAY RODENT BAIT


- Contains the single-feed bromadiolone, originally developed by Liphatech.
- Highly mould-resistant, yet extremely effective.
- 20g extruded block. Use less than other bromadiolone baits.
- Ideal all-year round maintenance bait.



## BEST USED:


- In commercial premises with ongoing rodent baiting programs to keep populations under control all year round
- Ideal for locations with high ambient humidity
- Can be used in a program to avoid resistance build-up in rodent populations.  
*Eg. alternate with a different bait containing a different active ingredient, such as Generation® Block.*
- Registered on Black rats, Norway rats and House mice, for use in and around buildings

### Comparative palatability of Maki Blocks against other baits (% consumption)



## HIGH MOULD-RESISTANCE

- Prolonged mould resistance for long-term efficacy



Tests run in a dairy factory with:  
15 °C and 87% relative humidity  
(ie: baiting 12 Jan 2007 and assessment 03 April 2007)

## WOULD YOU LIKE IT WRAPPED?

- Protect your bait from other pests (eg slugs, snails, insects, etc).
- Highly weather-resistant. Less melt, less mould.
- No direct contact with bait for end-user.
- Information printed on sachet makes on-site data recording easier.

## BEST USED:

- Mostly for commercial jobs, and for some residential situations (eg. subfloors, garden beds)
- Ideal when strong pressure from slugs/snails and/or insects (eg. cockroaches, etc.)
- Very appropriate for audited sites: each sachet displays some regulatory information making Maki Block Wrapped the product of choice for professional pest managers
- During summer months, when internal bait station temperatures get really high

## HOW TO USE:



Photo 1: bait preparation - Do not unwrap block.



Photo 2: to secure bait, fold plastic edge on the block and punch it neatly on the spike.



Photo 3: bait placement inside station - station stays cleaner till consumption.



Photo 4: rodent gnawing on bait.



Photo 5: counter visit.



Photo 6: dead rat => bait consumed.

# MAKI WRAPPED

WE'VE GOT YOU COVERED

New concept!



## BENEFITS:

- Plastic sachet increases longevity of baits protecting them from humidity and melting ==> less bait to replace
- Station stays cleaner, no mess at the bottom, other pests tend to not hibernate in the station

NOTE: Maki Block Wrapped will be performing best at times of low to medium rodent pressure. Replace with an unwrapped bait such as Maki Block or Generation Block when rodent pressure increases.